








# COMPATIBILITY GUIDE




-  **Fully assembled** Mass Spectrometer Bench.
-  **Fully movable** with lockable castor wheels.
-  Integrated vacuum pump(s) **noise enclosure**.
-  Fan **heat extraction** system.
-  Hoses & wires **grommets**.
-  Vacuum pump(s) **dampening vibration** system.
-  Overheating **temperature alarm**.




Delivered fully assembled



Fully movable



Integrated noise enclosure



Fan heat extraction system  
Hoses & wires grommets



Dampening vibration system



Overheating temperature alarm

# COMPATIBILITY GUIDE

TECHNOLOGY	BRUKER INSTRUMENTS	ionBench MODELS
MS only 	<b>EVOQ LC-TQ Qube, amaZon series</b>	BCH120NE58
	<b>EVOQ DART-TQ+ - short version</b>	BCH8097-DART-SET p/n <b>1907005</b>
	<b>EVOQ DART-TQ+ - large version</b>	BCH16597-DART-SET p/n <b>1907004</b>
	<b>EVOQ LC-TQ Elite</b>	BCH120NE78
	<b>compact, impact II VIP</b>	BCH120NECLH66_IM PACT-SET p/n <b>1907006</b>
	<b>timsTOF SCP, timsTOF Ultra</b>	BCH140-BRUKER-SET p/n <b>1887558</b>



**timsTOF SCP**  
BCH140-BRUKER-SET




**amaZon series**  
BCH120NE58



**EVOQ DART-TQ+**  
BCH8097-DART-SET

# COMPATIBILITY GUIDE

TECHNOLOGY	BRUKER INSTRUMENTS	ionBench MODELS
<p>LC/MS</p> 	<p><b>EVOQ LC-TQ Qube, amaZon series &amp; nanoElute &amp; Elute LC series</b></p>	<p>BCH190NE58</p>
	<p><b>EVOQ LC-TQ Elite &amp; nanoElute &amp; Elute LC series</b></p>	<p>BCH190NE78</p>
	<p><b>compact, impact II VIP &amp; nanoElute &amp; Elute LC series</b></p>	<p>BCH120NECLH66_ IMPACT-SET p/n <b>1907006</b></p>
	<p><b>timsTOF HT, timsTOF Pro &amp; nanoElute &amp; Elute LC series</b></p>	<p>BCHLC6575-BRASET p/n <b>1892593</b></p>



**nanoElute LC series &  
timsTOF Pro**  
BCHLC6575-BRASET




**Elute LC Series &  
amaZon series**  
BCH190NE58



**Elute LC series &  
EVOQ LC-TQ Elite**  
BCH190NE78

# COMPATIBILITY GUIDE

TECHNOLOGY	BRUKER INSTRUMENTS	ionBench MODELS
GC/MS 	<b>EVOQ GC-TQ Speed</b>	BCH120NE58
	<b>GC-APCI</b>	BCHLC_GC-APCI-SET p/n <b>1906944</b>



**EVOQ GC-TQ Speed**  
 BCH120NE58




**GC-APCI**  
 BCHLC\_GC-APCI-SET

# COMPATIBILITY GUIDE




- Fully assembled** Liquid Chromatography Bench.
- Fully movable** with lockable castor wheels.
- Electrically **height adjustable**.
- Chemical resistant** worksurface.
- Tubing management**.
- Minimize dead volume** between LC and MS.
- Safe access** to mobile phase bottles.




Delivered fully assembled



Fully movable



Weight capacity up to 250 kg



Tubing path system for waste management



Digital display for height adjustment



Chemically resistant and easy to clean.

# COMPATIBILITY GUIDE

TECHNOLOGY	BRUKER INSTRUMENTS	ionBench MODELS
LC only	<b>nanoElute &amp; Elute LC series</b>	BCHLC6575-BRASET p/n <b>1892593</b>
	<b>Elute LC series with Column oven</b>	BCHLC6765
	<b>Elute OLE</b>	BCHLC9075



**nanoElute LC series**  
BCHLC6575-BRASET



**Elute OLE**  
BCHLC9075



**Elute LC series with Column oven**  
BCHLC6765

## CONTACT US



### EUROPE

+33 (0) 9 70 44 65 57

[contact@ionbench.com](mailto:contact@ionbench.com)



### US

+1 984 219 8279

[contact\\_us@ionbench.com](mailto:contact_us@ionbench.com)