

COMPATIBILITY GUIDE




- Fully assembled** Mass Spectrometer Bench.
- Fully movable** with lockable castor wheels.
- Integrated vacuum pump(s) **noise enclosure**.
- Fan **heat extraction** system.
- Hoses & wires **grommets**.
- Vacuum pump(s) **dampening vibration** system.
- Overheating **temperature alarm**.




Delivered fully assembled



Fully movable



Integrated noise enclosure



Fan heat extraction system
Hoses & wires grommets




Dampening vibration system



Overheating temperature alarm

COMPATIBILITY GUIDE

TECHNOLOGY	PERKIN ELMER INSTRUMENTS	ionBench MODELS
LC/MS 	QSight SP50 & QSight LX50 & QSight220, QSight 420	BCH120NE78
	QSight SP50 & QSight220, QSight 420 & QS-Works Autosampler	BCH120NE78 GREY-COLOR




QSight SP50 & QSight LX50 & QSight 420
BCH120NE78



QSight SP50 & QSight 420 & QS-Works Autosampler
BCH120NE78 GREY-COLOR

COMPATIBILITY GUIDE

TECHNOLOGY	PERKIN ELMER INSTRUMENTS	ionBench MODELS
GC/MS 	GCMS 2400 System	BCH120NE58
	GCMS 2400 System & HS 2400 M Headspace Autosampler	BCH190NE58



GCMS 2400 System
BCH120NE58



**GCMS 2400 System & HS 2400 M
Headspace Autosampler**
BCH190NE58

COMPATIBILITY GUIDE

TECHNOLOGY	PERKIN ELMER INSTRUMENTS	ionBench MODELS
ICP-MS 	NexION 1000, NexION 2000B, NexION 5000 Multi-Quadrupole & noise enclosure for <i>vacuum pump</i> only	BCH190NE58
	NexION 1000, NexION 2000B, NexION 5000 Multi-Quadrupole & noise enclosure for <i>vacuum pump</i> only & <i>water chiller</i>	BCH190NE58NEPK
ICP-OES	Avio 220 Max, Avio 550 Max, Avio 560 Max & S23 Autosampler	BCH17097-AVIO



NexION 1000
BCH190NE58NEPK



Avio 220 Max & S23 Autosampler
BCH17097-AVIO




NexION 5000 Multi-Quadrupole
BCH190NE58

COMPATIBILITY GUIDE




- Fully assembled** Liquid Chromatography Bench.
- Fully movable** with lockable castor wheels.
- Electrically **height adjustable**.
- Chemical resistant** worksurface.
- Tubing management**.
- Minimize dead volume** between LC and MS.
- Safe access** to mobile phase bottles.




Delivered fully assembled



Fully movable



Weight capacity up to 250 kg



Tubing path system
for waste management



Digital display for height
adjustment



Chemically resistant and easy to
clean.

COMPATIBILITY GUIDE

TECHNOLOGY	PERKIN ELMER INSTRUMENTS	ionBench MODELS
LC only	LC 30, NexSAR, Flexar - 1 Stack	BCHLC4575
	LC 30, NexSAR - 1 Stack & Column Oven	BCHLC6765
	LC 30, NexSAR, Flexar - 2 Stacks	BCHLC9075



Flexar
BCHLC4575



NexSAR & Column Oven
BCHLC9075



LC30 UHPLC
BCHLC9075

CONTACT US



EUROPE

+33 (0) 9 70 44 65 57

contact@ionbench.com



US

+1 984 219 8279

contact_us@ionbench.com