








# COMPATIBILITY GUIDE




-  **Fully assembled** Mass Spectrometer Bench.
-  **Fully movable** with lockable castor wheels.
-  Integrated vacuum pump(s) **noise enclosure**.
-  Fan **heat extraction** system.
-  Hoses & wires **grommets**.
-  Vacuum pump(s) **dampening vibration** system.
-  Overheating **temperature alarm**.




Delivered fully assembled



Fully movable



Integrated noise enclosure



Fan heat extraction system  
Hoses & wires grommets



Dampening vibration system



Overheating temperature alarm

# COMPATIBILITY GUIDE

TECHNOLOGY	SCIEX INSTRUMENTS	ionBench MODELS
MS only	Triple Quad 3500, QTRAP 4500, Triple Quad 4500, Triple Quad 5500+ QTRAP Ready, X500B, X500R	BCH120NE58
	QTRAP 6500+, Triple Quad 6500+	BCH120NE78
	Triple Quad 7500 QTRAP Ready	BCH110100-SCIEX7500
	ZenoTOF 7600	BCH190NE58-LEFT
LC/MS	Triple Quad 3500, QTRAP 4500, Triple Quad 4500, Triple Quad 5500+ QTRAP Ready, X500B, X500R & ExionLC Series	BCH190NE58
	QTRAP 6500+, Triple Quad 6500+ & ExionLC Series	BCH190NE78
	ZenoTOF 7600 & ExionLC Series	BCH190NE58-LEFT & BCHLC6765
	Triple Quad 7500 QTRAP Ready & ExionLC Series	BCH19097-SCIEX7500



**Triple Quad 7500 QTRAP Ready  
&  
ExionLC Series**  
BCH19097-SCIEX7500



**ExionLC & ZenoTOF 7600**  
BCHLC6765 & BCH190NE58-LEFT



**QTRAP 6500+**  
BCH120NE78

# COMPATIBILITY GUIDE

TECHNOLOGY	SCIEX INSTRUMENTS	ionBench MODELS
MS & Dry pump	Triple Quad 7500 QTRAP Ready	BCH110100-SCIEX7500-DRY
	QTRAP 6500+, Triple Quad 6500+	BCH13097-SCIEX6500-DRY
LC/MS & Dry pump	QTRAP 6500+, Triple Quad 6500+ & ExionLC Series	BCH19097-SCIEX6500-DRY
	Triple Quad 7500 QTRAP Ready & ExionLC Series	BCH19097-SCIEX7500-DRY



**Triple Quad 6500 + & Dry pump**  
BCH13097-SCIEX6500-DRY



**ExionLC & Triple Quad 7500 QTRAP Ready & Dry pumps**  
BCH19097-SCIEX7500-DRY




**Triple Quad 7500 QTRAP**  
BCH13097-SCIEX7500-DRY

# COMPATIBILITY GUIDE



- Fully assembled** Liquid Chromatography Bench.
- Fully movable** with lockable castor wheels.
- Electrically **height adjustable**.
- Chemical resistant** worksurface.
- Tubing management**.
- Minimize dead volume** between LC and MS.
- Safe access** to mobile phase bottles.




Delivered fully assembled



Fully movable




Weight capacity up to 250 kg



Tubing path system for waste management



Digital display for height adjustment



Chemically resistant and easy to clean.

# COMPATIBILITY GUIDE

TECHNOLOGY	SCIEX INSTRUMENTS	ionBench MODELS
LC only	ExionLC, ExionLC 2.0, ExionLC 2.0+	BCHLC6765
	ExionLC AD	BCHLC9075



**ExionLC**  
BCHLC6765



**ExionLC AD**  
BCHLC9075

## CONTACT US



### EUROPE

+33 (0) 9 70 44 65 57

[contact@ionbench.com](mailto:contact@ionbench.com)



### US

+1 984 219 8279

[contact\\_us@ionbench.com](mailto:contact_us@ionbench.com)